



VIDAS® Easy SLM

Detection and confirmation of Salmonella in food products and environment

1 detection using VIDAS Easy SLM



2 confirmation with reference method



Salmonella range

42043	Buffered peptone water	6 x 225 ml
42629	Buffered peptone water	3 x 3 l
42118	Salmonella Xpress broth	20 tubes 10 ml
43621	SM ID2 - agar	20 plates
43563	XLD - agar	20 plates
43564	Hektoen - agar	100 plates
43111	Hektoen - agar	20 plates
20100	API 20 E	25 strips
32400	ID 32 E	25 strips
30702	VIDAS Salmonella (SLM)	60 tests

For further details, consult the product package insert.